

Transfection Introduction

What is Gene Delivery?

Gene delivery is the process of
introducing foreign genetic material (DNA or RNA) into cells
to modify or study their function.



1. What is Gene Delivery?

Types of genetic material insertion method



Transfection

Inserting DNA/RNA into
eukaryotic cells

Transformation

A **bacteria** take up DNA
and incorporate it into
their own genome

Transduction

Using a **virus** to transport
DNA/RNA into cells



1. What is Gene Delivery?


Step 1




Step 2




Step 3



Healthy cells with well-structured cell membrane



Temporary breakage of cell membrane
(e.g. electroporation)



DNA or RNA inserted into cells

Normally, cell does not absorb(allow) foreign substances without channels for specific substances with access. It only allows specific substances inside, which requires complex and strict steps.

Transfection is the key that makes this process a lot easier.



2. Purpose of Gene Delivery

Each cell produces proteins based on the genetic information.

This also means that we can manipulate cells by inserting the information we want to produce new protein or target protein.

Consideration 1

What would happen if we could choose the function of a cell?

Consideration 2

What if we could make a cell to only produce specific protein?
Or make a cell into something totally different?



2. Purpose of Gene Delivery

Generally to either manipulate cells to have specific genetic information or to study treatment



Purpose 1

For novel **drug development**
and disease treatment
research

Purpose 2

To **study gene expression**

Purpose 3

To **produce target protein**



2. Purpose of Gene Delivery



Novel drug development & research

- To study a way to fix the abnormal cells or to treat a genetic disease
- Widely used by inserting suicide gene inside cancer cells or inserting normal gene into a patient with genetic disorder
- E.g.) replacing old/not functioning parts into a new properly working parts



2. Purpose of Gene Delivery




Gene expression study

- To study how cells change/react or what cells produce when a certain type of gene is inserted.
- To study what the target gene's function is
- E.g.) pressing each button to see what it does



2. Purpose of Gene Delivery



Target protein production

- To manipulate cell to make massive amount of the target protein such as insulin or antibodies needed for drug
- E.g.) a factory producing proteins in bulk using a blueprint

Transfection

1. Broad Applicability in Eukaryotic Cells

- Transformation is limited to bacteria; transduction uses viral vectors with safety concerns

2. Safety

- Non-viral transfection reduces immunogenicity and insertional mutagenesis risk

3. Compatibility with Advanced Technologies

- Essential for CRISPR gene editing, mRNA therapeutics, and vaccine development


4. Cost-effectiveness and Ease of Use

- Simpler and more economical than viral vector production



1. Types of Transfection

Chemical Transfection




Chemical transfection uses

positively charged reagents to coat negatively charged nucleic acids,
facilitating their entry into cells by neutralizing charge repulsion.



1. Types of Transfection

Biological Transfection(Transduction)




Biological transfection(viral transduction) utilizes engineered viruses to deliver genetic material into host cells by leveraging the virus's natural ability to infect and insert genes into the genome.



1. Types of Transfection

Electroporation



Electroporation transfection delivers genetic material into cells by applying an **electrical pulse that temporarily permeabilizes the cell membrane.**



1. Types of Transfection

Types based on the purpose

	Chemical	Physical	Biological
Description	Using chemical substance to help penetrate through cell membrane	Using a device, a direct action given to a cell to make a temporary pores on the membrane	Using a virus, gene is transferred to a cell
Mechanism	1. Binds to gene and forms a complex 2. Enters the cell by fusion with the membrane or via endocytosis	1. Making temporary entrance on the membrane for genes to enter 2. Physically inserting the gene inside a cell	A virus is manipulated and used to transportation of gene
Example	Lipofection, calcium phosphate	Electroporation, microinjection	Adenovirus, lentivirus
Advantages	1. Simple 2. Cost-effective	Can be applied to various cell types	High efficiency and stability of gene transfer
Disadvantages	1. Efficiency may vary depending on the types of cell 2. Can be toxic to cells	1. Possible cell damage 2. Requires a device	1. Safety issue 2. Complex steps



1. Types of Transfection

Market Size

	Competitor	Market Size
Viral infection	Addgene, SBI, Sigma-Aldrich, Thermo-Fisher	Market share: 30% Revenue: 750 million USD
Lipofection	Thermo-Fisher (Lipofectamine), Polyplus, Mirus	Market share: 60% Revenue: 1.5 billion USD
Electroporation	Neon, Lonza, Bio-Rad, MaxCyte, BTX	Market share: 10% Revenue: 25 million USD


ExTransfection™

Electroporation system





History of ExTransfection




- **Proven Longevity** : Neon has been a trusted electroporation system for 16 years.
- **Global adoption** : Over **7,200 units** have been installed worldwide.
- **High consumable demand** : Annually, **15,000–20,000 consumables** are sold.
- **Research contributions** : Widely used in scientific studies and innovative projects.
- **Customer loyalty** : Maintains strong trust through quality and support.



ExTransfection_Technology

$$\text{Electric Field } E = \frac{V}{d}$$



ExTransfection™ Gold-Tip

Small electrode surface area with an **insulated capillary** ensures a consistent electric field due to the extended distance between the electrodes.


Bubble formation and electrolysis tends to be **minimized** with Gold-Tip. Thus, ExTransfection™ Gold-Tip **increases cell transfection efficiency**.

Consistent Electric Field

Cell damage is **minimized**

Transfection efficiency & cell viability are **enhanced**.

VS



Conventional Cuvette

The **short distance** between the electrodes and the **wide electrode surface area** result inconsistent electric field in cuvette.

Bubble formation and edge effects in cuvettes result in **uneven cell distribution**, leading to an **inconsistent electric field**.

Inconsistent Electric Field

Cell damage is **increased**


Transfection efficiency & cell viability are **lower**.



ExTransfection_Technology



ExTransfection™ Gold-Tip



Conventional Cuvette

Uniform electric field due to small electrode surface area

Electric Field

Uneven electric field due to large electrode surface area

Minimal heat generation

Heat Generation

Excess heat generated

Maintains **stable pH**

pH Stability

pH fluctuations due to electrode corrosion

Produces **fewer ions**

Ion Formation

Increased ion release

Gold-coated electrodes resist corrosion

Durability

Aluminum electrodes degrade over time

Low mechanical stress (high cell viability)

Shear Force

High mechanical stress




ExTransfection Kit Components

1. **Resuspension buffer R:** buffer commonly used to resuspend standard cell lines with DNA for electroporation.
2. **Resuspension buffer T:** buffer specifically designed for T cells and other sensitive cell types, optimized to maximize cell viability during electroporation.
3. **Electrolytic buffer E:** For use with 10uL ExTransfection tip
4. **Electrolytic buffer E2:** For use with 100uL ExTransfection tip






ExTransfection_How to use





ExTransfection_Demo Procedure

Cell line library



Electroporation !




1st Day





2nd Day

GFP, Flow cytometry




Analysis





ExTransfection_Demo Procedure





Key values and Differentiators

Customer feedback

"It was genuinely helpful. I can now perform tests I previously **avoided due to low transfection efficiency.**"

— Alina Mares

"It's **easy to use**, slightly more affordable than the previous system, and has great features like the database function. I highly recommend it."

— Dr. Lydia Wunderley

"**Very cost-effective. It always succeeds in transfection**, saving time and effort. It's incredibly efficient and reliable."

— Claire Demenis

"**It's much faster than reagent-based transfection systems.** Instead of taking over an hour, it completes transfection in just under 5 minutes."

— Mathlew Hayley



Neon™ NxT Electroporation System vs. Neon™ Transfection System

Specification	Neon NxT Electroporation System	Neon Transfection System
Pipette type	1-channel and 8-channel options	1-channel only
Electroporation Volume	10 µL or 100 µL	10 µL or 100 µL
Electroporation buffer volume*	2 mL	3 mL
Tip attachment	ClipTip technology	Friction
Electroporation pulses	1–10	1–10
Pulse duration	1–100 ms	1–100 ms
Pulse voltage	500–2,500 V	500–2,500 V
Arc detection	Yes	No
Cloud connectivity	Yes	No
Pulse generator dimensions**	9.5 x 7.6 x 9.9 in. (W x H x D) 11.9 lb (5.4 kg)	9.5 x 8.9 x 13.6 in. (W x H x D) 13.8 lb (6.25 kg)
Cable management feature***	Yes	No
Touch display	8-inch capacitive touchscreen	7-inch touchscreen
Electrical rating	100–240 VAC, 270 W	100–240 C, 150 W



- Pre-saved, frequently used cell-specific optimization protocols are stored.
- Tablet PC instead of touch screen