



Accessories for UVITEC Alliance Q9 Imaging System Series

TRANSILLUMINATORS

UV Transilluminators

Part Numbers	Description	Emission wavelength
2801 2602 1	Standard UV transilluminator, single wavelength, 26 cm x 21 cm	312 nm
2801 2606 1	Standard UV transilluminator, dual wavelength, 20 cm x 20 cm	312 nm/365 nm
2801 2600 1	Uvipure UV trans illuminator, single wavelength, 26 cm x 21 cm	312 nm
2801 2604 1	Uvipure UV transilluminator, dual wavelength, 20 cm x 20 cm	312 nm/365 nm

Blue Light Transilluminators

Part Numbers	Description	Emission wavelength
2801 2008 1	LED blue light transilluminator, 20 cm x 20 cm	470 nm

White Light Transilluminators

Part Numbers	Description
2801 2009 1	LED white light transilluminator, 20 cm x 20 cm

CONVERSION SCREENS

Part Numbers	Description
9120 2002 1	UV to white light conversion screen
9120 2003 1	UV to blue light conversion screen
9120 2004 1	Blue light to white light conversion screen

CHROMAPURE LIGHTING MODULES

Part Numbers	Description	Wavelength (nm)	Application
CLM-365	CHROMAPURE lighting module	UV 365 nm	TLC plates and Qdots
CLM-420	CHROMAPURE lighting module	Blue 420 nm	DyLight 405, CFP, eCFP, Cerulean, Atto 425 and Alexa 405
CLM-460	CHROMAPURE lighting module	Blue 460 nm	Cy2, GFP, eGFP, FITC, Alexa 488, Sypro, Bodipy, Atto 465 and DyLight 488
CLM-520	CHROMAPURE lighting module	Green 520 nm	Cy3, Cy3.5, RFP, Alexa 532 and 568, Krypton, Atto 550 and DyLight 549
CLM-650	CHROMAPURE lighting module	Red 650 nm	Cy5, Alexa 647, Atto 633, DyLight 633 and DyLight 649
CLM-690	CHROMAPURE lighting module	Red/NIR 690 nm	Cy5.5, Atto 680, IRdye 680, Alexa 680 and 700, and DyLight 680
CLM-740	CHROMAPURE lighting module	NIR 740 nm	Cy7, Alexa 750 and DyLight 750
CLM-780	CHROMAPURE lighting module	IR 780 nm	IRdye 800, Alexa 790 and DyLight 800

NARROW BAND-PASS FILTERS

Part Numbers	Description	Peak wavelength (nm)	Bandpass Range (nm)
1710 4123 1	Narrow band-pass filter (F-520)	520 nm	± 25 nm from peak wavelength
1710 4180 1	Narrow band-pass filter (F-535)	535 nm	± 25 nm from peak wavelength
1710 4160 1	Narrow band-pass filter (F-565)	565 nm	± 25 nm from peak wavelength
1710 4181 1	Narrow band-pass filter (F-595)	595 nm	± 25 nm from peak wavelength
1710 4161 1	Narrow band-pass filter (F-655)	655 nm	± 25 nm from peak wavelength
1710 4182 1	Narrow band-pass filter (F-695)	695 nm	± 25 nm from peak wavelength
1710 4162 1	Narrow band-pass filter (F-710)	710 nm	± 25 nm from peak wavelength
1710 4184 1	Narrow band-pass filter (F-740)	740 nm	± 25 nm from peak wavelength
1710 4186 1	Narrow band-pass filter (F-760)	760 nm	± 25 nm from peak wavelength
1710 4183 1	Narrow band-pass filter (F-820)	820 nm	± 25 nm from peak wavelength
1710 4185 1	Narrow band-pass filter (F-850)	850 nm	± 25 nm from peak wavelength

SAMPLE TRAYS

Part Numbers	Description
9111 1158 1	Chemiluminescence sample plate for western blot imaging
9111 1162 1	Fluorescence sample plate for transparent samples (e.g. gels)